

Test Description

NRAS mutation test is an in vitro diagnostic test for the qualitative detection of mutations in codons 12, 13, 59, 61, 117 and 146 of NRAS oncogene.

Patient Demographic

Name: Dinesh Singh Yadav
Sex: Male
Date of Birth/Age: 51 Years
Disease: Colorectal Cancer

Clinician

Clinician Name: Dr Sandeep Batra
Medical Facility: Max Hospital
Pathologist: Not Provided

Specimen

Site: Rectal Sigmoid Growth
Sample Type: FFPE block B 5852/19
Date of Collection: 11-05-2019
Date of Booking: 11-05-2019

NRAS Mutation Analysis

Result

No Mutation Detected in NRAS

GENOMIC FINDINGS

No mutation detected

INTERPRETATION

No Mutation detected in NRAS codons 12, 13, 59, 61, 117 and 146.

METHODOLOGY

The NRAS Mutation Test, performed on the Biocartis Idylla™ system, is an *in vitro* diagnostic test for the qualitative detection of 18 mutations (G12C, G12S, G12D, G12A, G12V, G13D, G13V, G13R, A59T, Q61H/Q61H, Q61K/R/L, K117N/K117N and A146T/V) in codons 12, 13, 59, 61, 117 and 146 of the NRAS gene. Formalin-fixed paraffin-embedded (FFPE) human cancer tissue is lysed liberate DNA for subsequent real-time PCR amplification using allele specific primers. Two sample processing controls (SPC) are amplified simultaneously i.e. (1) a conserved region of the NRAS gene (referred to as NRAS-Total) and (2) a conserved region of the BRAF gene. The presence of a mutant genotype is determined by calculating the difference between the NRAS Sample Processing Control Cq and the Cq obtained for the NRAS mutant signal(s).

The analytic sensitivity of this assay has been determined at < or = 5%

REFERENCES

1. Allegra et al. J Clin Oncol (2016) 34:179-85
2. Boleij et al. BMC Cancer (2016) 16:825.



May 13, 2019

Dr Gulshan Yadav, MD, Consultant Pathology

Date