


RAL  
TECNICA PARA EL LABORATORIO

BIO RAD

05/02/2024

N=624

SEROPROTEINS



PRINTED IN U.S.A.

TOT. PROTEINS: 7.8

A/G = 1.2

FRACTION	%	G/DL	% NORM.
ALBUMIN	54.0	4.2	52 - 68
ALPHA1	2.3	0.2	2 - 5
ALPHA2	9.8	0.8	8 - 12
BETA	16.6	1.3	10 - 16
GAMMA	17.4	1.4	12 - 18

BIO RAD